

# Hospital Furnitures



# Hospital Furniture

## HOSPITAL BED STEAD GENERAL SWE-110210

- Fixed height bed
- Overall size 1950 mm X 910 mm X 600 mm
- Sheet top
- Legs with rubber feet
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- Supplied in knock down construction

### Optional

SWE-110210-C

With Castors (100mm dia.); Set of 4, Two with brakes.



## HOSPITAL BED STEAD GENERAL SWE-110214

- Fixed height
- Overall size 1950 mm X 910 mm X 600 mm
- Rigid mesh platform
- Legs with rubber feet
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- Supplied in knock down construction

### Optional

SWE-110214-C

With Castors (100mm dia.); Set of 4, Two with brakes



## HOSPITAL BED STEAD GENERAL SWE-110110

- Overall size 1970 mm X 910 mm X 600 mm
- Perforated sheet top
- Removable head & foot bows of 31 mm dia. with six vertical support on head side & 4 vertical support on foot side.
- Legs with rubber feet
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- Knock down construction

### Optional

SWE-110110-C

With Castors (100mm dia.); set of 4, two with brakes



# Hospital Furniture

## HOSPITAL BED STEAD-Dx SWE-110310

- Overall size 1970 mm X 910 mm X 600 mm
- Perforated Sheet top
- Removable Stainless Steel tubular head & foot bows with laminated panels
- Legs mounted on Castors (100mm dia); two with brakes
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- Supplied in knock down construction



## HOSPITAL BED STEAD TWO SECTION, MANUAL SWE-111210

- Backrest adjustable, Manual
- Overall size 1950 mm X 910 mm X 600 mm
- Two section Rigid mesh top
- Legs with rubber feet
- Pre-treated epoxy powder coated
- Provision for Mosquito net poles.
- Knock down construction

### Optional

SWE-111210-C

With Castors (100mm dia.); set of 4, Two with brakes

## HOSPITAL BED STEAD TWO SECTION, MANUAL SWE-111120

- Backrest adjustable, Manual
- Overall size 1970 mm X 910 mm X 600 mm
- Two section perforated sheet top.
- Removable head & foot bows (31mm dia.) on both sides with six vertical support on head side, 1060mm height and 4 vertical supports leg side. 820mm height
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- knock down construction

### Optional

SWE-111120-C

With Castors (100mm dia.); set of 4, Two with brakes



# Hospital Furniture

## HOSPITAL BED STEAD TWO SECTION, MANUAL SWE-111310

- Backrest adjustable, Manual
- Overall size 1970 mm X 910 mm X 600 mm
- Perforated sheet top.
- Removable Stainless Steel tubular head & foot bows with laminated panels
- Legs mounted on Castors (100mm dia.); two with brakes
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- Supplied in knock down construction



## HOSPITAL BED STEAD TWO SECTION, MECHANICAL SWE-111230

- Backrest adjustable through rotating crank
- Overall size 1950 mm X 910 mm X 600 mm
- Two section Rigid mesh top
- 100mm four swivelling castors; 2 with brakes
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- Supplied in knock down construction

## HOSPITAL BED STEAD TWO SECTION, MECHANICAL SWE-111240

- Backrest adjustable through crank
- Overall size 1970 mm X 910 mm X 600 mm
- Perforated sheet top
- Stainless Steel I.V. rod with four locations
- Legs with rubber feet
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- Knock down construction

### Optional

SWE-111240-C

With Castors (100mm dia.); set of 4, two with brakes



# Hospital Furniture

## HOSPITAL BED STEAD TWO SECTION, MECHANICAL SWE-111150

- Backrest Adjustable through Rotating Crank
- Overall size 1970 mm X 910 mm X 600 mm
- Perforated Sheet top
- Four corner rubbers buffers
- Stainless Steel I.V. rod with four locations
- Mounted on 100 mm swivelling castors, two with brakes
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- Supplied in knock down construction



## HOSPITAL BED STEAD TWO SECTION, MECHANICAL SWE-111340

- Backrest adjustable through crank
- Overall size 1970 mm X 910 mm X 600 mm
- Perforated sheet top
- Four corner rubbers buffers
- Stainless Steel I.V. rod with four locations
- Removable Stainless Steel tubular head & foot bows with laminated panels
- 100mm four swivelling castors; 2 with breaks
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- Supplied in knock down construction

## BED STEAD FOUR SECTION, FOWLER POSITION SWE-112210

- Overall size 1970 mm X 910 mm X 600 mm
- Adjustable backrest & leg section by means of two independent cranks
- Perforated sheet top
- Head bows 1060mm (H) and foot bows 820mm (H) with three vertical tubular supports
- 100mm four swivelling castors; 2 with breaks
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- Knock down construction

### Optional

SWE-112210-RM

Bed Stead four section rigid mesh platform



# Hospital Furniture

## HOSPITAL BED STEAD FOUR SECTION, FOWLER POSITION

SWE-112112

- Backrest and knee break position adjustable through crank
- Overall size: 1970 mm X 910 mm X 600 mm
- Perforated sheet top
- Four corner rubbers buffers
- Stainless Steel I.V. rod with four locations
- Removable tubular head & foot bows
- 100 mm of swivelling castors ; 2 with brakes
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- knock down construction

### Optional

SWE-112112-RM

With rigid mesh platform



## BED STEAD FOUR SECTION, FOWLER POSITION

SWE-112300

- Bed Stead, Four sections
- Overall size 1970 mm X 910 mm X 600 mm
- Perforated sheet top.
- Adjustable backrest & leg section by means of two independent cranks
- Four corner rubber buffer
- Stainless Steel I.V. rod with four location
- Removable Stainless Steel tubular head & foot bows with laminated panels
- 100mm four swivelling castors 2 with brakes
- Pre - treated epoxy powder coated
- Provision for Mosquito net poles.
- Knock down construction

### Optional

SWE-112300 - RM


With rigid mesh platform



## BED FOR HOME CARE

SWE-111500

- Overall Size: 2100 x 1100 x 450 mm
- Base, frame and leg made of Mild Steel Epoxy Powder coated
- Head-foot boards & side panels made of quality wood.
- Backrest adjustable by cabled remote control
- Castors: 125mm, two with brakes



# Hospital Furniture



## INTENSIVE CARE BED

SWE-114333

### Overall Size:

Head & Foot Bows

- 2110mm (L) x 910mm (W) x 500-750 mm (H)
- Backrest: 0-80°
- Knee-Rest: 0-40°
- Trendelenburg & Reverse: 0-12°

### Frame Work:

- Frame Made of Rectangular Steel Tube with Anti-Bacterial Epoxy Powder Coating

### Platform:

- Four Section platform made of Zinc Coated, Mild Steel, Epoxy Powder Coated
- Stainless Steel Wooden Laminated Tubular Head & Foot Bows, easily removable

### Side Rails:

- Dropdown side rails, stainless steel, on both sides

### Castors:

- Four Swivelling Castors 125mm, two with brakes

### Operations:

- Backrest, Head Rest, Knee Rest, Trendelenburg/Reverse and through Four Separate cranks
- Provided with corner Buffers

### Provisions:

- I.V. Pole (Telescopic) at four locations
- Two locations hooks for urine bag holder

### SWE-113330

Intensive Care Bed-Fixed Height

- Overall Size: 2250 x 1015 x 720mm
- Four Section Platform with Trendelenburg & Reverse Trendelenburg Position

# Hospital Furniture

## INTENSIVE CARE BED, ELECTRIC

SWE-114444



CE

### Overall Size:

- | 2110 mm (L) x 910 mm (W) x 450-750 mm (H)
- | Backrest: 0-80°
- | Knee-Rest: 0-40°
- | Trendelenburg & Reverse: 0-12°
- | Static Load Capacity: 240Kg.

### Frame Work:

- | Frame Made of Rectangular Steel Tube with Anti-Bacterial Epoxy Powder Coating

### Platform:

- | Four Section platform made of Zinc Coated Mild Steel, Latest parallel beam type surface
- | Epoxy Powder Coated

### Head & Foot Boards

- | Strong ABS engineering plastic easily removable

### Side Rails

- | Four Half Length tuck-away side rails, on both side

### Castors

- | Four Swivelling Castors 125 mm, two with Brakes

### Operations:

- | Backrest, Knee Rest, Trendelenburg/ Reverse and Height Adjustment through cabled remote control.

### Provisions:

- | I.V. Pole (Telescopic) at four locations
- | Two locations hooks for urine bag holder

### SWE-113440

#### Intensive Care Bed-Fixed Height

- | Overall Size: 2250 x 1015 x 720mm
- Backrest, Knee Rest and Trendelenburg/Reverse through cabled remote control and manual override.

# Hospital Furniture

## BED ACCESSORIES



ABS PLASTIC HEAD AND FOOT PANELS

ACC-101011

COLLAPSIBLE ALUMINIUM SIDE RAILS

ACC-101012



CENTRAL BRAKING CASTORS

ACC-101013

TUCK AWAY (SPLIT) SIDE RAILS SMALL

ACC-101014



TUCK AWAY (SPLIT) SIDE RAILS BIG

ACC-101015

COLLAPSIBLE TELESCOPIC SIDE RAILS

ACC-101016



# Hospital Furniture

**BALKAN FRAME**

**ACC-101017**



**BABY CRIB WITH ATTACHEMENT**

**ACC-101018**



**BED ELEVATOR**

**ACC-101019**



**LIFTING POLE**

**ACC-101020**



**MOSQUITO NET POLES**

**ACC-101021**



**I.V. ROD**

**ACC-101022**



# Hospital Furniture

## BED SIDE LOCKER

SWE-115110

- Overall size 400 (L) X 400 (W) X 820 mm (H)
- Stainless Steel top
- Machine pressed Mild steel box with door having ventilation lowers IS 4266:1967
- Self locking knob type latched, mounted on rubber shoes
- Pre-treated Epoxy powder coated



CE



CE

## BED SIDE LOCKER

SWE-115220

- Overall size 400 (L) X 400 (W) X 820 mm (H)
- Three side covered having one drawer & cupboard mounted on 5 cm swivelling castors
- Stainless steel top with raised edges on three



CE

## BED SIDE LOCKER

SWE-115225

- Overall size 400 X 400 X 820 mm (H)
- Bed side Locker with one drawer and full cupboard
- Stainless steel top
- Pre-treated epoxy powder coated

# Hospital Furniture



## BED SIDE LOCKER

**SWE-115200**

- Laminated top
- Overall size 400 (L) X 400 (W) X 820 mm (H)
- Three side machine pressed steel sheet having one drawer and cupboard mounted on 2.5cm swivelling castors at the rear
- Front of drawer, cupboard & top made of laminated wood
- Pre-treated Epoxy powder coated

## BED SIDE LOCKER

**SWE-115600**

- Dimensions: 480 x 480 x 760 mm
- Stylish design, Light weight, fully made of ABS plastic
- Lifting bars and towel clips on both side of the cabinet
- Absolutely inert to stains and chemicals. Very easy to clean
- Drawer with easy movement and cupboard with removable shelves with bottle holder.
- Sliding tray for multi purposes and with glass holder



## BEDSIDE LOCKER WITH OVERBED TABLE

**SWE-116100**

- Overall size : 650 mm (L) x 400 mm (W) x 850 mm (H)
- Three side machine pressed steel sheet
- Having one drawer & cupboard mounted on four 2.5cm swivelling castors
- Laminated Front, drawer & cabinet
- Laminated over bed table; height adjustable from 90cm to 120cm.

# Hospital Furniture

## OVER BED TABLE, TWO SECTION TOP

SWE-116200

- Overall size : 810 mm (L) x 355 mm (W) x 800-1050 mm (H)
- Height adjustable through rotating crank
- Tubular steel frame work
- Two section laminated top, for easy reading
- Four 50mm swivelling castors.
- Pre-treated Epoxy Powder coated
- Knock down construction

### Optional:

SWE-116200-PH, Height adjustable through pneumatic pump



## OVERHEAD TABLE, HEIGHT ADJUSTABLE

SWE-116311

- Overall size : 810 mm (L) x 355 mm (W) x 800-1050 mm (H)
- Height adjustable
- Tubular steel frame work
- Wooden Laminated/ABS top
- Four 50mm swivelling castors.
- Pre-treated Epoxy Powder coated
- Knock down construction

### Optional:

SWE-116311-PH, Height adjustable through pneumatic pump with ABS top (pic. B)



## OVER BED TABLE, FIXED HEIGHT

SWE-116400

- Overall size : 1060 mm (L) x 360 mm (W) x 1000 mm (H)
- Table top size : 850mm (L) x 370mm (W)
- Fixed Height
- Tubular steel frame work
- Laminated top
- Four 50mm swivelling castors.
- Pre-treated Epoxy Powder coated
- Knock down construction



# Hospital Furniture

## REVOLVING STOOL

SWE-142150

- Stainless Steel top
- Height adjustable from 450mm to 680mm, through threaded screws
- Four legged base made of 25mm steel tube mounted on rubber shoes
- Stainless Steel ring for footrest.
- Pre-treated Epoxy powder coated frame work

**SWE -142159**

Revolving Stool

- Complete Stainless Steel Construction



CE



CE

(A)



CE

(B)

## REVOLVING STOOL

SWE-142100

- Cushioned top
- Height adjustable from 450 mm to 680 mm, through threaded screws
- Four legged made of 25mm steel tube mounted on rubber shoes
- Stainless Steel ring for footrest.
- Pre-treated Epoxy powder coated

**Optional:**

New design patient stool with back support and castors (pic. B)

## REVOLVING STOOL

SWE-142200

- Height adjustable through pneumatic pump
- Five spoke Plastic, unbreakable base, with castors
- Stainless steel ring for footrest

**Optional:**

With S.S. Top (pic. B)



(A)



(B)

CE

# Hospital Furniture



**REVOLVING STOOL**

**SWE-143100**

- Plastic Seat & Back
- Adjustable height 500 to 700 mm, through pneumatic pump
- Five legged base on castors



**REVOLVING STOOL**

**SWE-143200**

- Cushioned seat & back with adjustable armrest
- Height Adjustment through hydraulic pump
- Five legged base on castors,
- Stainless steel foot support for foot rest



CE



**ATTENDANT COUCH CUM CHAIR**

**SWE-119300**

- Overall size:  
Chair position: 770 x 680 x 1125 mm  
Bed position: 2050 x 680 x 550 mm
- Stainless steel tubular frame work mounted
- Easily convertible to chair
- Castors to facilitate couch position
- High quality durable mattress in three sections
- Compact construction to save space

# Hospital Furniture



**BED SIDE SCREEN, TWO FOLDS SWE-144200**

- 2 Section Tubular Steel Frame, Epoxy Powder Coated
- Washable Plastic Curtains
- Mounted on Castors

**BED SIDE SCREEN, THREE FOLDS SWE-144300**

- Tubular steel frame, Epoxy Powder Coated
- Mounted on Castors
- Washable Plastic Curtains
- Knock Down Construction



**BEDSIDE SCREEN, FOUR FOLDS SWE-144400**

- Tubular steel frame, Epoxy Powder Coated
- Mounted on Castors
- Washable Plastic Curtains
- Knock Down Construction

# Hospital Furniture

## FOOT STEP SINGLE

**SWE-146100**

- Overall size : 510 mm X 300 mm X 230 mm (H)
- Top with rubber matt/aluminum tread flats.
- Surrounded by angle frame.
- Tubular frame work.
- Pre-treated Epoxy powder Coated.

**SWE-146199**

Complete Stainless Steel



## FOOT STEP DOUBLE

**SWE-146200**

- Overall size : 510 mm X 460 mm X 460 mm
- Two Steps
- Top with rubber matt/aluminum tread flats.
- Surrounded by angle frame.
- Tubular frame work.
- Pre-treated Epoxy powder Coated.

**SWE-146299**

Complete Stainless Steel

## MAYO INSTRUMENT TROLLEY

**SWE-131100**

- Overall size : 560 mm X 400 mm X 760 mm (H)
- Height adjustable
- Stainless steel tray size 560 mm (L) X 400mm (W)
- Tubular frame work of mild steel: epoxy Powder coated; mounted on 50 mm swivelling castors.
- Knock down construction.

**SWE-131199**

Complete Stainless Steel



# Hospital Furniture



## MAYO INSTRUMENT TROLLEY, DOUBLE ENDED

**SWE-131200**

- Overall size : 1100 mm X 600 mm X 1170 - 1850 mm (H)
- Height adjustable
- Stainless steel tray size 1100 mm (L) X 600mm (W)
- Tubular frame work mild steel epoxy powder coated
- Powder coated mounted on 75 mm swivelling castors.
- Knock down construction.

**SWE-131299**

Complete Stainless Steel

## INSTRUMENT TROLLEY WITH RAILING

- Stainless Steel shelves with three sides guard rails
- Epoxy powder Coated tubular frame work.
- Mounted on 750mm swivelling castors.

**SWE-133648**

600mm (L) X 450mm (W) X 800 mm (H)

**SWE-133768**

750 mm (L) x 600mm (W) X 800mm (H)

**SWE-133968**

900 mm (L) x 600mm (W) X 800mm (H)



## ALL STAINLESS STEEL CONSTRUCTION

**SWE-134648**

600mm (L) X 450mm (W) X 800 mm (H)

**SWE-134768**

750 mm (L) x 600mm (W) X 800mm (H)

**SWE-134968**

900 mm (L) x 600mm (W) X 800mm (H)



## INSTRUMENT TROLLEY

- Two Stainless Steel shelves with upstands on four sides
- One Push Handle
- Complete stainless steel
- Mounted on 750mm swivelling castors.

**SWE-135648**

600mm (L) X 450mm (W) X 800 mm (H)

**SWE-135768**

750 mm (L) x 600mm (W) X 800mm (H)

**SWE-135968**

900 mm (L) x 600mm (W) X 800mm (H)


# Hospital Furniture

## OPTIONAL ACCESSORIES FOR INSTRUMENT TROLLEYS



Complete Trolley

CE



Basin Holder with  
Stainless Steel Basin

CE



Removable Waste Bin,  
Stainless Steel

CE



Stainless Steel Bottle Holder  
for 6 Bottles

CE



One Drawer

CE



Two Drawers

CE

# Hospital Furniture



## DRESSING CARRIAGE

**SWE-135600**

- Overall size : 750 x 450 x 800 mm
- With Two Stainless Steel Shelves with guard rails on three sides
- Side Railing With Bracket For Bucket and Bowl
- Mild steel tubular frame work epoxy powder coated
- Mounted on 750 mm Swivelling castors

**SWE-135699**

Complete Stainless Steel


## DRESSING TROLLEY WITH DRAWERS

**SWE-135722**

- Two drawers
- Overall size : 610 X 460 mm
- S.S Shelve with three side raised edges
- Mild steel tubular frame work, epoxy powder coated
- Mounted on 750 mm Swivelling castors

**SWE-135799**

Complete Stainless Steel



## MEDICINE TROLLEY

**SWE-190161**

- Overall Size: 40 x 60 x 110 cm
- Stainless Steel Construction
- Upper Section with a large open space
- Six Drawers and a storage compartment with shelf and locking door
- Mounted on castors

CE

# Hospital Furniture

## BABY CRADLE

SWE-167065

- Swinging Type
- Size: Stand 770 mm (L) X 600 mm (W) X 880 mm (H)
- Cradle 750mm (L) X 400mm (W) X 770mm (H)
- Tubular steel stand on four rubber stumps
- Pretreated & Epoxy powder Coated



## PAEDIATRIC BED WITH DROP SIDE RAILS

SWE-167060

- Paediatric bed with drop side rails
- Overall size: 1480mm (L) X 840mm (W) X 600mm (H)
- M.S. Sheet top, both bows of 1060 mm (H) with one horizontal support and six vertical supports
- Provision for Mosquito net poles
- Legs fitted with 75 mm castors
- Full length drop side rails
- Pretreated & Epoxy powder Coated

## BABY CRIB

SWE-167067

- Trolley manufactured from stainless steel tube
- Mounted on four castors 60mm.
- Crib manufactured from Plexiglass.
- Trendelenburg & Reverse Trendelenburg position.
- Crib Size: 820 x 460 x 260mm



# Hospital Furniture



CE

## SWADDLING TABLE

SWE-167069

- Overall Size : 1150 (L) X 650 (W) X 900 mm (H)
- Cushioned top covered with washable material
- Frame made of Mild Steel, epoxy Powder coated
- Laminated lower shelf
- Knock down construction

## SWADDLING TABLE

SWE-167068

- Padded swaddling table mounted on Metallic cupboard



CE



CE

## GYNECOLOGICAL COUCH

SWE-153110

- Overall size : 1800 mm x 550 mm x 820 mm
- Made of rectangular Mild steel tube
- Adjustable back rest & leg section
- Top is upholstered & covered with washable material, flame retardant
- Pair of lithotomy crutches
- Stainless Steel Bowl
- Knock down construction

SWE-153199

Complete Stainless Steel

# Hospital Furniture

## BABY CRADLE

SWE-167065

- Swinging Type
- Size: Stand 770 mm (L) X 600 mm (W) X 880 mm (H)
- Cradle 750mm (L) X 400mm (W) X 770mm (H)
- Tubular steel stand on four rubber stumps
- Pretreated & Epoxy powder Coated



## PAEDIATRIC BED WITH DROP SIDE RAILS

SWE-167060

- Paediatric bed with drop side rails
- Overall size: 1480mm (L) X 840mm (W) X 600mm (H)
- M.S. Sheet top, both bows of 1060 mm (H) with one horizontal support and six vertical supports
- Provision for Mosquito net poles
- Legs fitted with 75 mm castors
- Full length drop side rails
- Pretreated & Epoxy powder Coated



## BABY CRIB

SWE-167067

- Trolley manufactured from stainless steel tube
- Mounted on four castors 60mm.
- Crib manufactured from Plexiglass.
- Trendelenburg & Reverse Trendelenburg position.
- Crib Size: 820 x 460 x 260mm

CE

# Hospital Furniture

## DELIVERY TABLE

SWE-153150

- Overall size : 1800 mm x 550 mm x 820 mm
- Made of Mild steel rectangular steel tube
- Three sections with independent controls and inclinable seat section
- Pair of lithotomy crutches
- Stainless Steel Bowl
- Knock down construction

SWE-153199

Complete Stainless Steel



## DELIVERY BED

SWE-152110

- Overall size : 1800 x 550 mm x 820 mm (H)
- Made of Mild Steel square tube
- Two section, leg section moveable on castors; with two brakes
- Leg section can slide under the main frame
- Cushioned upholstered top covered with washable flame retardant material.
- Pair of lithotomy crutches
- Stainless Steel Bowl
- Knock down construction.

SWE-152199

Complete Stainless Steel



## DELIVERY BED

SWE-153120

- Overall size : 1800 x 550 mm x 820 mm (H)
- Made of Mild Steel square tube
- Three section, leg section movable on castors; with two brakes
- Leg section can slide under the main frame
- Cushioned upholstered top covered with washable flame retardant material.
- Backrest Adjustable
- Pair of lithotomy crutches
- Stainless Steel Bowl
- Knock down construction.

SWE-153199

Complete Stainless Steel

# Hospital Furniture



**OBSTETRIC LABOUR TABLE**

**SWE-153200**

- Overall size : 1870 x 600 x 800 mm
- Mild Steel tubular frame work, Stainless Steel top in three section
- Backrest adjustable on ratchet
- Middle section with perineal cut out
- Leg section can slide under main frame
- Gas spring operated Trendelenburg & Reverse Trendelenburg Position
- Stainless Steel sides railings
- Lithotomy leg holders with straps, one set of hand grip
- Provision for Stainless Steel bowl
- Pre-treated Epoxy powder coated

**Optional**

Cushioned lithotomy crutches

SWE-153299: Complete Stainless Steel

**OBSTETRIC LABOUR TABLE**

**SWE-153300**

- Overall size : 1870 x 600 x 800 mm
- Mild Steel tubular frame work
- Upholstered top in three section
- Backrest adjustable on ratchet
- Leg section can slide under main frame
- Trendelenburg Position
- Lithotomy leg holders with straps
- One set of hand grip
- Provision for Stainless Steel bowl
- Pre-treated Epoxy powder coated



**OBSTETRIC LABOUR TABLE**

**SWE-153400**

- Overall size : 1870 x 600 x 800 mm
- Mild Steel tubular frame work
- Upholstered top in three section
- Height Adjustable Through Foot Operated Hydraulic Pump
- Backrest adjustable on ratchet
- Leg section can slide under main frame
- Trendelenburg Position
- Padded lithotomy crutches
- Provision for Stainless Steel bowl
- Pre-treated Epoxy powder coated



# Hospital Furniture

## OBSTETRIC LABOUR TABLE

**SWE-152500**

- Overall size : 1870 x 600 x 800mm
- Mild Steel Tubular Framework with upholstered top in three sections.
- Back rest, trendelenburg and height adjustable through cabled remote control
- Leg section can slide under main frame
- Padded lithotomy crutches
- Provision for Stainless Steel bowl
- Pre-treated Epoxy powder coated



## STRETCHER TROLLEY

**SWE-121100**

- Overall size 2140 x 510 x 820mm (H)
- Tubular frame work mounted on four swivelling castors 150mm dia.
- Removable stretcher top made of mild steel sheet
- Stainless Steel tubular handles
- Provision for I.V. rod
- Pre-treated Epoxy powder coated
- Knock down construction

### Optional

- Cylinder cage
- Side rails
- One end raised on ratchet
- Cushioned upholstered top, Covered with washable material

SWE-121199  
complete Stainless Steel

## EXAMINATION TABLE TWO FOLD

**SWE-126100**

- Overall size: 1830 x 510 x 820 mm (H)
- Two fold, backrest adjustable on ratchet
- Tubular frame work, fitted with rubber shoes
- Pre-treated Epoxy Powder Coated
- Knock down construction
- Set of Mattress



# Hospital Furniture



**EXAMINATION TABLE, CUSHIONED TOP**

**SWE-126200**

- Overall size: 1830 x 510 x 820 mm (H)
- Two fold, backrest adjustable on ratchet
- Cushioned top
- Mild steel tubular frame work, fitted with rubber shoes

**Optional:**

Lithotomy strips  
Paper roll holder

**EXAMINATION BED**

**SWE-126300**

- Size: 1830 x 510 x 820 mm
- Backrest 0-40°, adjustable by gas spring
- Two sectioned cushioned top
- Three drawers on telescopic channels and cupboard
- BP machine holder near backrest
- Mild steel pre-treated epoxy powder coated finish
- Sturdy construction for longevity
- Washproof 2 section fixed mattress of high quality



**PLASTER TABLE**

**SWE-187189**

- Complete Stainless Steel Construction
- Top section with 6 moveable panels
- Dimension 180X65X80 cm
- Container (for fluid and waste) adjustable at any panel.

# Hospital Furniture

## EMERGENCY RECOVERY TROLLEY

SWE-121300

- Overall size 2150 X 710, mattress size 1900 mm X 550 mm, two section top
- Height adjustable by foot operated hydraulic pump from 660 to 910 mm
- X-ray permeable top with X-ray cassette holder
- Telescope IV rod
- Stainless Steel side rails
- Backrest adjustable on ratchet
- Trendelenburg & reverse trendelenburg positions through gas spring
- Mounted on four swivelling castors of 150 mm dia. two with breaks

Optional

Oxygen cylinder holder, tray.



## EMERGENCY RECOVERY TROLLEY/ STRETCHER

SWE-121500

- Dimensions Overall: 2120 x 800 X 610-930 mm
- Stretcher Size: 1830 (L) x 710 mm (W)
- Static Load Capacity 180 Kgs
- Castors: Single lever, Central and unique direction lock
- Four section top made of ABS Engineering Plastic with adjustable backrest and fowler position
- Radio translucent top with X-ray cassette holder
- Collapsible aluminium side rails
- Foot operated Trendelenburg / Reverse adjustment
- Hydraulic assisted foot operated height adjustment
- 4 smooth swivelling castors of 200 mm with Central Braking System (CBS)
- Rectilinear treadle for straight motion and Steering handles on the head side of the patient for enabling single attendant manoeuvrability
- Space for stowing Oxygen cylinder at the bottom of the trolley
- Buffers at all 4 corners
- Provision for attaching IV stand on all 4 corners
- Mattress retainers on both sides of the patient
- Rubberized push bar for ease of pushing and absorption of impact
- Sturdy design, Superior finish for life time of service
- ABS covered easy to clean base as per international design

# Hospital Furniture



CE

## OXYGEN CYLINDER TROLLEY

SWE-187152

- Suitable for 10 litre size gas cylinder
- Mild steel tubular frame work fitted with two wheels of 100mm dia.
- Pre-treated epoxy powder coated.

## CRASH CART

SWE-197175

- Overall size : 940 mm x 500 mm x1540 mm (H)
- Mild steel tubular framework
- Six coloured removable bins & two polystyrene
- Lockable storage units with three drawers each
- Four swiveling castors of 125mm dia., two with brakes
- Corner rubber buffers
- Oxygen cylinder holder
- Electric examination lamp
- Cardiac massage board
- Stainless Steel I.V. rod
- Wooden laminated shelves
- Knock down construction

SWE-197199

Complete Stainless Steel



CE



CE



CE

## WASH BASIN: DOUBLE

SWE-187137

- Overall Size : Two basins each of 350mm Dia. X 850mm (H)
- Mild steel with Epoxy powder coated tubular Frame work mounted on 55mm swiveling castors.
- Stainless steel 350mm basin

SWE-187199

Complete stainless steel

# Hospital Furniture

## SOLID LINEN BIN

**SWE-187138**

- Overall Size : 910 x 510 mm (Dia.)
- Mild steel tubular framework
- Mounted on 75 mm dia. castors with detachable canvas bag.

**SWE-187199**

complete Stainless Steel



## KICK BUCKET

**SWE-187239**

- Overall Size : 300 dia x 300 (H)
- Mild steel Epoxy powder Coated frame mounted on 75mm swiveling castors
- Stainless steel container

**SWE-187299**

Complete Stainless Steel

## SWAB RACK

**SWE-187305**

- Overall Size: 137 (H) x 76 (W) x 46 (D) cms
- Stainless Steel Tubular Framework, mounted on 10 cms Dia Castor.
- Eight Rows; Six Hooks in each Row.
- Stainless Steel tray fitted at Bottom



# Hospital Furniture



## I.V. STAND

### SWE-187950

- Stainless Steel I.V. Rod with two locations
- Mild Steel Epoxy Powder Coated column
- Moulded base with five branches on castors

### SWE-187959

- Stainless Steel I.V. Rod and column

### SWE-187950

- Stainless Steel I.V. Rod with four locations
- Mild Steel Epoxy Powder Coated column
- Moulded base with five branches on castors

### SWE-187953

- Stainless Steel I.V. Rod with two locations
- Mild Steel Epoxy Powder Coated column
- Moulded base with three branches, on castors

# Hospital Furniture


## INSTRUMENT & MEDICINE CABINET

SWE-197185

- Made of Mild Steel Sheet
- Lockable glass doors
- Adjustable Stainless Steel/ Glass shelves
- Overall Size : 1800mm X 800mm X 400mm

SWE-197851

- Overall size : 1450mm X 600mm X 350mm



## INSTRUMENT & MEDICINE CABINET

SWE-193300

- Overall size : 1800 x 800 x 400mm
- Mild Steel Construction
- Divided in two section, upper section with two sliding glass doors & with two SS/glass shelves
- Lower section with two metal doors with locking
- Removable two SS/glass shelves
- Pretreated Epoxy powder coated



## INSTRUMENT & MEDICINE CABINET

SWE-193333

- Overall size : 1800 x 800 x 400 mm
- Mild Steel Construction
- Divided in two sections, Both with two sliding glass doors & with two SS/glass shelves
- Pretreated Epoxy powder coated

# Hospital Furniture



CE

## WHEEL CHAIR

SWE-171001

### Non Folding

- 425mm Cushioned seat & back, Fixed armrests
- Front with 125mm castors and rear 650mm solid wheels, with breaks
- Pr-treated Epoxy powder coated

SWE - 171099

Stainless Steel Construction

## WHEEL CHAIR

SWE-172002

- Folding
- Removable cushioned seat & adjustable cushioned back
- Removable armrest
- Adjustable leg rest with cushion pad support
- Adjustable Foot rest
- Front with 150mm solid castors & rear 650mm solid wheels with breaks
- Pr-treated Epoxy powder coated

SWI - 172099

Stainless Steel Construction



CE



CE

## BLOOD DONOR CHAIR

SWE-120300

- Overall size : 1800 mm (L) x 600 mm (W)
- Cushioned padded three sections top covered with washable material
- Backrest and leg rest adjustable simultaneously
- Stainless steel tubular frame work

# Hospital Furniture

## STRETCHER FOR AMBULANCE

SWE-128200

- Size : 1900 x 560 x 900 mm (in standing mode)  
1900 x 560 x 280 mm (in folded mode)
- Made of high quality alloyed Aluminium
- Top surface made from foamed cushion
- Foldable legs controlled by levers on foot end.
- Legs lockage by chute mechanism
- Weight : 40 Kgs.
- 125 mm Casters for effortless movement
- Load bearing capacity : 140 Kgs.



## COLLAPSIBLE CHAIR CUM TROLLEY STRETCHER

SWE-128300

- Size: 190 x 56 x 80-85 cm
- Weight: 30 kg.
- Adjustable back rest 0°-90°
- Loading Capacity 160-180 kg.
- Automatic loading stretcher with capability to convert into wheel chair.
- Anodized aluminum construction
- Completely overturnable Side protections

## FOLDING WHEEL CHAIR CUM STRETCHER

SWE-128700

- Dimensions: 690 x 510 x 910 mm
- Folding Dimensions: 180 x 510 x 910 mm
- Aluminium frame work
- Simple to operate
- Light weight 8.5 kg. Nett.
- Loading capacity: 159 kg
- Manufactured from durable material, easy to clean & wash




# Hospital Furniture

## SCOOP STRETCHER

SWE-128100

- Portable
- Design: Open:  $2000 \pm 10$  mm  
Closed:  $1650 \pm 10$  mm  
Head Position:  $420 \pm 10$  mm  
Foot Position:  $305 \pm 05$  mm
- Divide locking device
- Aluminium alloy for tubular framework.
- 3 Velcro straps 4 inches in which to immobilize head, chest and lower limbs.



## STRETCHER FOLDING DOUBLE FOLD

SWE-128500

- Tube: Aluminium 12 Gauge Thickness
- Locking & Folding Clamps: Stainless Steel
- Cloth: Canvas of good quality
- Load Bearing: 120 kg.
- Carrying Bag: Canvas of good quality
- Dimensions (Unfolded size): 2130 (L) x 560 (W) x 150 (H) mm



## STRETCHER FOLDING

SWE-128400

- Tube: Mild Steel
- Locking & Folding Clamps: Stainless Steel
- Cloth: Canvas of good quality
- Load Bearing: 120 kg.
- Dimensions: 2130 L x 560 W x 150 H mm
- Easily foldable for easy stowage



## STRETCHER FOLDING FOUR FOLD

SWE-128600

- Tube: Aluminium 12 Gauge Thickness
- Locking & Folding Clamps: Stainless Steel
- Cloth: Canvas of good quality
- Load Bearing: 150 kg.
- Carrying Bag: Canvas of good quality
- Dimensions: Unfolded size: 220 x 53 x 15 cm



**Works:**

1701-04, HSIDC, Industrial Estate, Rai, Distt: Sonapat, Haryana (INDIA)  
Tel: +91-130-2367562 (10 lines) Fax: +91-130-2367561

**Corporate & Sales Office:**

B-69, Industrial Area, G.T. Karnal Road, Delhi-110033 (INDIA)  
Tel: +91-11-4708 3333 (8 lines)  
Fax: +91-11-4709 3333  
E-mail: [surgicoin@vsnl.com](mailto:surgicoin@vsnl.com), Web: [www.surgicoin.com](http://www.surgicoin.com)

